

# AlloSeq Tx Genotyping Using DNA Extracted from Saliva

Supaneda N. Kolanski<sup>1</sup>, Peter B. Dallas<sup>1</sup>, Hayley M. Hogan<sup>1</sup>, Jessica A. Etccl<sup>1</sup>, Christopher N. Newbound<sup>1</sup>, Lindsey T. Madden<sup>1</sup>, Samantha L. Timler<sup>1</sup>, Emma JD. Carr<sup>1</sup>, Damian M. Goodridge<sup>1</sup>, and David C. Sayer<sup>1</sup>

<sup>1</sup>CareDx Pty Inc, Fremantle, Australia

## INTRODUCTION

The AlloSeq<sup>®</sup> Tx Hybrid Capture assay is a flexible and comprehensive methodology for the typing of human leukocyte antigen (HLA) genes, and other important non-HLA transplant-associated genes, at high resolution. Although the AlloSeq<sup>®</sup> Tx assay together with AlloSeq<sup>®</sup> Assign<sup>™</sup> software enables researchers to robustly genotype genomic DNA samples extracted from whole blood, the compatibility of the system for genotyping DNA from saliva samples has not been assessed. In order to assess the performance of this sample type, 24 genomic DNA samples extracted from saliva were tested using the AlloSeq<sup>®</sup> Tx assay.

## METHODS & MATERIALS

24 individuals provided saliva samples using Oragene DNA OG-600 (saliva) and ORAcollect DNA OCR-100 (sponge) kits (DNA genotek). Genomic DNA was subsequently extracted with a QIAamp DNA Mini Kit (Qiagen). Library construction and Enrichment were undertaken using the AlloSeq<sup>®</sup> Tx system. During this experiment, in-development Hybridisation Capture probes targeting 19 loci were utilised including HLA-A, -B, -C, -DRB1, -DQB1, -DQA1, -DPB1, -DPA1, -DRB3/4/5, -E, -F, -G, -H, MICA, MICB and two additional genes. Enriched libraries were then loaded onto a micro flowcell on the Illumina MiSeq System for Next Generation Sequencing. Finally, AlloSeq<sup>®</sup> Assign<sup>™</sup> software was used for data analysis and allele assignments of the 19 target loci.




Figure 1. AlloSeq<sup>®</sup> Tx assay workflow

## RESULTS

The average concentration of genomic DNA extracted from the Oragene saliva and sponge samples was 15.32 ng/ $\mu$ L and 6.78 ng/ $\mu$ L respectively, above the minimum sample input required for the AlloSeq<sup>®</sup> Tx assay. AlloSeq<sup>®</sup> Tx libraries were generated from the genomic DNA extracted from both saliva and sponge samples with an average yield of 24.93 ng/ $\mu$ L and an average fragment size of 788 bp. The average enrichment pool yield was 14.05 ng/ $\mu$ L with an average fragment size of 728 bp. The resulting sequence data exhibited 94% of reads with an average read quality of Q30 or higher, an average read depth of 140x, and average allele balance of 42%. High resolution genotyping results were achieved with 4 field typing for Class I loci, 3 field typing results for Class II and 2 fields typing for MICA/MICB for all samples tested. All genotyping results were 100% concordant with QTYPE<sup>®</sup> genotyping results.




Figure 2. Library concentration yield and distribution achieved from the Oragene saliva samples compared to the yield from Oragene sponge samples. All samples were above the minimum input requirements for the AlloSeq<sup>®</sup> Tx assay.




Figure 3. TapeStation trace illustrating the fragment peak size and distribution of AlloSeq<sup>®</sup> Tx libraries produced from the Oragene saliva and sponge. The average peak height is 788bp.

## RESULTS (cont'd)

Field	HC161L121SAL22*	HC161L121SAL23	HC161L121SAL24	HC161L121SAL25	HC161L121SAL27	HC161L121SAL28	HC20177L136A01SM02	HC20177L136A02SM05	HC20177L136A03SM06	HC20177L136A04SM07	HC20177L136A05SM09	HC20177L136A06SM10	HC20177L136B01SM12	HC20177L136B02SM13	HC20177L136B03SM14
IMGT/A	01:01:01:01	07:02:01:01	02:02:02:--	32:01:01:01	27:02:01:--	07:02:01:--	X	04:01:01G	01:02:01	05:02:01	15:01:01	03:01:01G	01:02:01	06:02:01	16:02:01
IMGT/B	02:01:01:01	38:01:01:01	06:02:01:01	01:03:01	04:01:01G	01:03:01G	03:03:01	01:02:01G	01:03:01	03:03:01	03:03:01	03:03:01	03:03:01	03:03:01	03:03:01
IMGT/C	26:01:01:01	57:01:01:01	12:03:01:01	X	06:01:01G	02:01:01	06:03:01	07:01:01	02:01:01	03:01:01	07:01:01	04:01:01	05:01:01	05:01:01	05:01:01
IMGT/DPA1	02:01:01:01	15:01:01	03:03:01:01	X	04:01:01	01:02:01	05:03:01	13:02:01	14:54:01	16:02:01	03:01:01	04:01:01	05:01:01	05:01:01	05:01:01
IMGT/DPB1	11:01:01:01	38:02:01	03:02:02:05	01:03:01	04:01:01	01:02:01G	02:01:01	16:02:01	03:01:01	03:01:01	03:01:01	03:01:01	03:01:01	03:01:01	03:01:01
IMGT/DQA1	33:03:01:01	58:01:01:03	07:02:01:01	02:02:02	05:01:01	05:01:01G	05:02:01	03:01:01	03:01:01	03:01:01	03:01:01	03:01:01	03:01:01	03:01:01	03:01:01
IMGT/DOB1	03:01:01:01	07:02:01:01	02:02:02:01	01:03:01	04:01:01	01:02:01	03:01:01	15:01:01	06:02:01	08:03:01	08:03:01	08:03:01	08:03:01	08:03:01	08:03:01
IMGT/DRB1	02:01:01:01	15:01:01:01	03:03:01:01	X	04:01:01G	01:03:01	03:01:01	01:01:01	04:01:01	03:01:01	03:01:01	03:01:01	03:01:01	03:01:01	03:01:01
IMGT/DRB3	24:02:01:01	44:02:01:01	05:01:01:02	X	04:02:01G	03:03:01	05:01:01	04:01:01	03:02:01	05:01:01	04:01:01	04:01:01	04:01:01	04:01:01	04:01:01
IMGT/DRB4	29:02:01:01	44:03:01:--	16:01:01:01	X	04:01:01	01:02:01	02:02:01	15:01:01	06:02:01	07:01:01	07:01:01	07:01:01	07:01:01	07:01:01	07:01:01
IMGT/DRB5	02:01:01:01	45:01:01:03	16:01:01:01	02:01:08	01:01:01	01:01:02	05:01:01	01:02:01	01:02:01	01:02:01	01:02:01	01:02:01	01:02:01	01:02:01	01:02:01
IMGT/DRB6	02:01:01:01	45:01:01:03	X	04:02	06:01:01	01:02:01	06:02:01	15:01:01	06:02:01	06:02:01	06:02:01	06:02:01	06:02:01	06:02:01	06:02:01
IMGT/DRB7	01:01:01:01	07:02:01:01	07:01:01:01	01:03:01	01:01:01	01:01:01	02:01:01	15:01:01	02:01:01	02:01:01	02:01:01	02:01:01	02:01:01	02:01:01	02:01:01
IMGT/DRB8	02:01:01:01	08:01:01:01	07:02:01:03	02:01:02	04:01:01	05:01:01	06:02:01	03:01:01	06:02:01	06:02:01	06:02:01	06:02:01	06:02:01	06:02:01	06:02:01
IMGT/DRB9	02:01:01:01	45:01:01:01	03:03:01:01	01:03:01	04:01:01G	01:04:01	02:02:01	14:54:01	02:02:01	02:02:01	02:02:01	02:02:01	02:02:01	02:02:01	02:02:01
IMGT/DRB10	26:01:01:01	55:01:01:01	06:02:01:03	X	06:01:01G	02:01:01	05:03:01	07:01:01	05:03:01	05:03:01	05:03:01	05:03:01	05:03:01	05:03:01	05:03:01
IMGT/DRB11	11:01:01:01	27:05:02:--	01:02:01:01	01:03:01	04:01:01	01:01:01	01:01:01	01:01:01	01:01:01	01:01:01	01:01:01	01:01:01	01:01:01	01:01:01	01:01:01
IMGT/DRB12	24:02:01:01	39:06:02:01	07:02:01:01	X	04:01:01	01:02:01	06:02:01	15:01:01	01:02:01	01:02:01	01:02:01	01:02:01	01:02:01	01:02:01	01:02:01
IMGT/DRB13	01:01:01:01	08:01:01:01	02:02:02:01	01:03:01	01:01:01	03:03:01	02:01:01	03:01:01	02:01:01	02:01:01	02:01:01	02:01:01	02:01:01	02:01:01	02:01:01
IMGT/DRB14	31:01:02:01	27:05:02:09	07:01:01:01	02:01:02	04:01:01	05:01:01	03:02:01	04:01:01	03:02:01	03:02:01	03:02:01	03:02:01	03:02:01	03:02:01	03:02:01
IMGT/DRB15	02:01:01:01	07:02:01:01	05:01:01	01:03:01	03:01:01G	03:01:01	03:01:01	04:01:01	03:01:01	03:01:01	03:01:01	03:01:01	03:01:01	03:01:01	03:01:01
IMGT/DRB16	03:01:01:01	44:02:01:01	07:02:01	X	04:01:01G	03:03:01	03:02:01	04:07:01	03:02:01	03:02:01	03:02:01	03:02:01	03:02:01	03:02:01	03:02:01
IMGT/DRB17	01:01:01:01	08:01:01:01	07:01:01:01	01:03:01	04:01:01	05:01:01	02:01:01	03:01:01G	05:01:01	05:01:01	05:01:01	05:01:01	05:01:01	05:01:01	05:01:01
IMGT/DRB18	25:01:01:01	58:01:01:01	07:18:01:01	X	X	05:05:01	03:01:01	13:03:01G	03:01:01	03:01:01	03:01:01	03:01:01	03:01:01	03:01:01	03:01:01
IMGT/DRB19	02:01:01:01	07:02:01:01	07:02:01:03	01:03:01	03:01:01	01:02:01	03:02:01	15:01:01	01:02:01	01:02:01	01:02:01	01:02:01	01:02:01	01:02:01	01:02:01
IMGT/DRB20	07:02:01:01	X	07:02:01:03	02:01:01	13:01:01	02:01:01	06:02:01	07:01:01	06:02:01	06:02:01	06:02:01	06:02:01	06:02:01	06:02:01	06:02:01

Figure 4. Screenshot of the summary view from the AlloSeq<sup>®</sup> Assign<sup>™</sup> software displaying 4 fields genotypes for Class I loci, and 3 fields for Class II loci. A number of novel alleles with variants in non-coding and/or UTR regions of class I loci were identified, which were reported to 3 fields.




Figure 5. Comparison of the %Q30 score observed for samples extracted from Oragene saliva and sponge samples compared with HLA Reference Standards (IHWG) controls processed via AlloSeq<sup>®</sup> Tx system. All samples' Q-score pass illumina's run specifications.




Figure 6. Comparison of the average coverage (or read depth) observed for samples extracted from Oragene saliva and sponge samples compared with HLA Reference Standards (IHWG) controls processed via AlloSeq<sup>®</sup> Tx system. A single outlier sample was observed with an average coverage of less than 100x, which still genotyped successfully.




Figure 7. Comparison of the average allele balance observed for samples extracted from Oragene saliva and sponge samples compared with HLA Reference Standards (IHWG) controls processed via AlloSeq<sup>®</sup> Tx system.

## CONCLUSIONS

This study has demonstrated successful high-resolution genotyping and high-quality data achieved from genomic DNA extracted from Oragene saliva and sponge samples with the AlloSeq<sup>®</sup> Tx assay. Performance with this sample type enables wider research application of the AlloSeq<sup>®</sup> Tx assay.

## ACKNOWLEDGEMENTS

Fremantle R&D team  
Connie Xu and CareDx AB Stockholm Research team